

Product Datasheet

Recombinant Human Oncostatin-M/OSM Protein (orb1551435)

Description	Recombinant Human Oncostatin-M/OSM Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.01 EU/ug by LAL method.
Target	OncostatinM
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1052, Recombinant Human Oncostatin-M/OSM Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala26-Arg221) of human Oncostatin-M/OSM (Accession #NP_065391.1).
Expiration Date	6 months from date of receipt.