

Product Datasheet

Recombinant Human EphA2 Protein (orb1551436)

Description	Recombinant Human EphA2 Protein
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	EphA2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1051, Recombinant Human EphA2 Protein is produced by insect Cell-baculovirus expression system. The target protein is expressed with sequence (Asp596-Ile976) of human EphA2 (Accession #NP_004422.2).
Expiration Date	6 months from date of receipt.