

## Product Datasheet

# Recombinant Human MAPK1 Protein, GST His Tag (orb1551438)

<b>Description</b>	Recombinant Human MAPK1 Protein, GST His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact us for more information.
<b>Target</b>	MAPK1
<b>Note</b>	For research use only
<b>Purity</b>	> 70% by SDS-PAGE.
<b>Source</b>	KMP1049, Recombinant Human MAPK1 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Ser360) of human MAPK1 (Accession #NP_002736.3) fused with a GST His Tag tag at the N- terminus.
<b>Expiration Date</b>	6 months from date of receipt.