

Product Datasheet

Recombinant Human GMF-beta Protein (orb1551441)

Description	Recombinant Human GMF-beta Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	GMFbeta
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1046, Recombinant Human GMF-beta Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser2-His Tag142) of human GMF-beta (Accession #NP_004115.1).
Expiration Date	6 months from date of receipt.