

Product Datasheet

Recombinant Human RBBP-4 Protein (orb1551442)

Description	Recombinant Human RBBP-4 Protein
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	RBBP4
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1045, Recombinant Human RBBP-4 Protein is produced by insect Cell-baculovirus expression system. The target protein is expressed with sequence (Met1-Ser425) of human RBBP-4 (Accession #NP_005601.1).
Expiration Date	6 months from date of receipt.