

Product Datasheet

Recombinant Human IL-36 alpha Protein, His Tag (orb1551447)

Description	Recombinant Human IL-36 alpha Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	IL36 alpha
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1040, Recombinant Human IL-36 alpha Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Phe158) of human IL-36 alpha (Accession #NP_055255.1) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.