

Product Datasheet

Recombinant Human Thioredoxin Protein, His Tag (orb1551448)

Description	Recombinant Human Thioredoxin Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Thioredoxin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1039, Recombinant Human Thioredoxin Protein is produced by E. coli expression system. The target protein is expressed with sequence (Val2-Val105) of human Thioredoxin (Accession #NP_003320.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.