

Product Datasheet

Recombinant Human SULT1A1 Protein, His Tag (orb1551450)

Description	Recombinant Human SULT1A1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.01 EU/ug by LAL method.
Target	SULT1A1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1037, Recombinant Human SULT1A1 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Leu295) of human SULT1A1 (Accession #NP_001046.2) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.