

Product Datasheet

Recombinant Human TNFSF13B Protein (orb1551451)

Description	Recombinant Human TNFSF13B Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	TNFSF13B
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1036, Recombinant Human TNFSF13B Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala134-Leu285) of human TNFSF13B (Accession #NP_006564.1).
Expiration Date	6 months from date of receipt.