

Product Datasheet

Recombinant Human VEGF 121 Protein (orb1551460)

Description	Recombinant Human VEGF 121 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	VEGF 121
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1027, Recombinant Human VEGF 121 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala207-Arg327 (Lys321Asn)) of human VEGF 121 (Accession #NP_001020541.2).
Expiration Date	6 months from date of receipt.