

Product Datasheet

Recombinant Human UBE2D3 Protein, His Tag (orb1551461)

Description	Recombinant Human UBE2D3 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	UBE2D3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1026, Recombinant Human UBE2D3 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Met147) of human UBE2D3 (Accession #NP_003331.1) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.