

## Product Datasheet

# Recombinant Human Serpin F1 Protein (orb1551462)

<b>Description</b>	Recombinant Human Serpin F1 Protein
<b>Species/Host</b>	Virus
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	Serpin F1
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1025, Recombinant Human Serpin F1 Protein is produced by insect Cell-baculovirus expression system. The target protein is expressed with sequence (Gln20-Pro418) of human Serpin F1 (Accession #NP_002606.3).
<b>Expiration Date</b>	6 months from date of receipt.