

## Product Datasheet

### Recombinant Human IL-18 Protein (orb1551463)

<b>Description</b>	Recombinant Human IL-18 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	IL18
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1024, Recombinant Human IL-18 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Tyr37-Asp193) of human IL-18 (Accession #NP_001553.1).
<b>Expiration Date</b>	6 months from date of receipt.