

Product Datasheet

Recombinant Human SLC3A2 Protein (orb1551464)

Description	Recombinant Human SLC3A2 Protein
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	SLC3A2
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1023, Recombinant Human SLC3A2 Protein is produced by insect Cell- baculovirus expression system. The target protein is expressed with sequence (Arg206-Ala630) of human SLC3A2 (Accession #NP_002385.3).
Expiration Date	6 months from date of receipt.