

## Product Datasheet

### Recombinant Human PPIA Protein (orb1551465)

<b>Description</b>	Recombinant Human PPIA Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	PPIA
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1022, Recombinant Human PPIA Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Glu165) of human PPIA (Accession #NP_066953.1).
<b>Expiration Date</b>	6 months from date of receipt.