

Product Datasheet

Recombinant Human S100A12 Protein, His Tag (orb1551469)

Description	Recombinant Human S100A12 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	S100A12
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1018, Recombinant Human S100A12 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Glu92) of human S100A12 (Accession #NP_005612.1) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.