

Product Datasheet

Recombinant human CD40/TNFRSF5 Protein, Fc & His Tag (orb1551474)

Description	Recombinant human CD40/TNFRSF5 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD40
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1013, Recombinant human CD40/TNFRSF5 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Glu21-Arg193) of human CD40/TNFRSF5 (Accession #NP_001241.1) fused with a Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.