

Product Datasheet

Recombinant Human LILRB2/CD85d/ILT4 Protein, His Tag (orb1551475)

Description	Recombinant Human LILRB2/CD85d/ILT4 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	LILRB2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1012, Recombinant Human LILRB2/CD85d/ILT4 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln 22 - Val 461) of human LILRB2/CD85d/ILT4 (Accession #AAB87662) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.