

Product Datasheet

Recombinant Human BTLA Protein, Fc & His Tag (orb1551480)

Description	Recombinant Human BTLA Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BTLA
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1007, Recombinant Human BTLA Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Lys31-Thr134) of human BTLA (Accession #NP_001078826.1.) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.