

Product Datasheet

Recombinant Human TNF-alpha Protein, His Tag (orb1551485)

Description	Recombinant Human TNF-alpha Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	TNF alpha
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1002, Recombinant Human TNF-alpha Protein is produced by E. coli expression system. The target protein is expressed with sequence (Val77-Leu233) of human TNF-alpha (Accession #NP_000585.2) fused with an initial Met at the N-terminus and a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.