

Product Datasheet

Human PRKAA1 Protein, His & GST Tag (orb1551490)

Description	Human PRKAA1 Protein, His & GST Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug determined by the LAL method
Target	PRKAA1
Note	For research use only
Purity	> 90%, determined by SDS-PAGE
Source	This protein contains the human PRKAA1 (Q13131-1) (Met 1-Gln 559) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus, constructed the plasmid 1; This protein contains the human PRKAB1 (Q9Y478) (Met 1-Ile 270), constructed the plasmid 2; This protein contains the human PRKAG1 (P54619) (Met 1-Pro 331), constructed the plasmid 3 . The three plasmids were co-expressed and the heterotrimer was purified and expressed from Baculovirus-Insect Cells.
Expiration Date	6 months from date of receipt.