

## Product Datasheet

### Human PFK1 Protein, His & GST Tag (orb1551494)

<b>Description</b>	Human PFK1 Protein, His & GST Tag
<b>Species/Host</b>	Virus
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug determined by the LAL method
<b>Target</b>	PFK1
<b>Note</b>	For research use only
<b>Purity</b>	> 90%, determined by SDS-PAGE
<b>Source</b>	This protein contains the human PFKM (P08237-1) (Thr2-Val780) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus and expressed from Baculovirus-Insect Cells.
<b>Expiration Date</b>	6 months from date of receipt.