

Product Datasheet

Human PPM1A Protein, His Tag (orb1551495)

| | |
|------------------------|---|
| Description | Human PPM1A Protein, His Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Target | PPM1A |
| Note | For research use only |
| Purity | > 85%, determined by SDS-PAGE |
| Source | This protein contains the human PPM1A isoform alpha-1 (P35813-1) (Met 1-Trp 382) was fused with a polyhistidine tag at the C-terminus and expressed from E. coli. |
| Expiration Date | 6 months from date of receipt. |