

## Product Datasheet

### Human MIF Protein, His Tag (orb1551497)

<b>Description</b>	Human MIF Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug determined by the LAL method
<b>Target</b>	MIF
<b>Note</b>	For research use only
<b>Purity</b>	> 92%, determined by SDS-PAGE
<b>Source</b>	This protein contains the human MIF (P14174) (Pro 2-Ala 115) was fused with a polyhistidine tag at the C-terminus and a signal peptide at the N-terminus and expressed from HEK293 Cells.
<b>Expiration Date</b>	6 months from date of receipt.