

Product Datasheet

Influenza A H5N1 (A/VietNam/1203/2004) HA Protein (HA1 Subunit), His & Fc Tag (orb1551500)

Description	Influenza A H5N1 (A/VietNam/1203/2004) HA Protein (HA1 Subunit), His & Fc Tag
Species/Host	Human
Endotoxins	< 1.0 EU/ug determined by the LAL method
Target	H5N1
Note	For research use only
Purity	> 90%, determined by SDS-PAGE and SEC-HPLC
Source	This protein contains an N-terminal segment (Met 1-Arg 342) of Influenza A Virus H5N1 Hemagglutinin (AAW807171), termed as HA1, was fused with the C-terminal His-tagged Fc region of mouse IgG1 at the C-terminus and expressed from HEK293 Cells.
Expiration Date	6 months from date of receipt.