

Product Datasheet

Recombinant Human ACE2 Protein, His-tag (orb1551503)

Description	Recombinant Human ACE2 Protein, His-tag
Species/Host	Chinese Hamster
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	ACE2
Note	For research use only
Purity	≥95% as determined by SDS-PAGE
Source	A DNA sequence encoding the human ACE2 (NP_068576.1) (Gln18-Ser740) was expressed with a polyhistidine tag at the C-terminus.
Expiration Date	6 months from date of receipt.