

Product Datasheet

Recombinant Human TNFRSF10C C-fc Protein, Fc & His Tag (orb1551515)

Description	Recombinant Human TNFRSF10C C-fc Protein, Fc & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	TNFRSF10C
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF10C(Ala26-Ala221) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	O14798
Expiration Date	6 months from date of receipt.