

Product Datasheet

Recombinant Human TNFa/TNFSF2/TNF Protein, No Tag (orb1551523)

Description	Recombinant Human TNFa/TNFSF2/TNF Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TNF
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNF(Val77-Leu233) was fused without Tag and expressed in E. coli.
Uniprot ID	P01375
Expiration Date	6 months from date of receipt.