

Product Datasheet

Recombinant Human TNF RII/TNFRSF1B Protein, His Tag (orb1551524)

Description	Recombinant Human TNF RII/TNFRSF1B Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TNFRSF1B
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF1B(Pro24-Thr206) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P20333
Expiration Date	6 months from date of receipt.