

Product Datasheet

Recombinant Human TNF RI/TNFRSF1A Protein, His Tag (orb1551525)

Description	Recombinant Human TNF RI/TNFRSF1A Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TNFRSF1A
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF1A(Ile22-Thr211) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P19438
Expiration Date	6 months from date of receipt.