

## Product Datasheet

# Recombinant Human TMX2/TXNDC14 Protein, His Tag (orb1551526)

<b>Description</b>	Recombinant Human TMX2/TXNDC14 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	TMX2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human TMX2(Met125-Lys296) was fused with the N-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q9Y320">Q9Y320</a>
<b>Expiration Date</b>	6 months from date of receipt.