

Product Datasheet

Recombinant Mouse TL1A/Tnfsf15 Protein, No Tag (orb1551529)

Description	Recombinant Mouse TL1A/Tnfsf15 Protein, No Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Tnfsf15
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Tnf15(Ile76-Leu252) was fused without Tag and expressed in E. coli.
Uniprot ID	Q5UBV8
Expiration Date	6 months from date of receipt.