

Product Datasheet

Recombinant Human TIM-3/HAVCR2 C-fc Protein, Fc & His Tag (orb1551539)

Description	Recombinant Human TIM-3/HAVCR2 C-fc Protein, Fc & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	HAVCR2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TIM3(Ser22-Arg200) was fused with the C-terminal Fc and His Tag and expressed in Mammalian cells.
Uniprot ID	Q8TDQ0
Expiration Date	6 months from date of receipt.