

Product Datasheet

Recombinant Mouse Thrombopoietin/Tpo/Thpo Protein, His Tag (orb1551548)

Description	Recombinant Mouse Thrombopoietin/Tpo/Thpo Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Thpo
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Tpo(Ser22Thr356) was fused with the N-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P40226
Expiration Date	6 months from date of receipt.