

## Product Datasheet

# Recombinant Human TGF- $\alpha$ Protein, No Tag (orb1551555)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human TGF- $\alpha$ Protein, No Tag   |
| <b>Species/Host</b>    | E. coli   |
| <b>Reactivity</b>      | Human   |
| <b>Endotoxins</b>      | Please contact our customer service for more information.   |
| <b>Target</b>          | TGFA  |
| <b>Note</b>            | For research use only   |
| <b>Purity</b>          | $\geq 90\%$ as determined by SDS-PAGE   |
| <b>Source</b>          | This protein contains the human TGF- $\alpha$ (Val41-Val90) was fused without Tag and expressed in E. coli. |
| <b>Uniprot ID</b>      | <b>P01135</b>   |
| <b>Expiration Date</b> | 6 months from date of receipt.  |