

Product Datasheet

Recombinant Mouse TGF-beta RII/Tgfbr2 Protein, Fc Tag (orb1551557)

Description	Recombinant Mouse TGF-beta RII/Tgfbr2 Protein, Fc Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Tgfbr2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse TGFR2(Ile24-Asp159) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	Q62312
Expiration Date	6 months from date of receipt.