

Product Datasheet

Recombinant Human TAF4 Protein, His Tag (orb1551569)

Description	Recombinant Human TAF4 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TAF4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TAF4(Ser35-Arg140) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q96LR4
Expiration Date	6 months from date of receipt.