

Product Datasheet

Recombinant Human SWSAP1 Protein, His Tag (orb1551578)

Description	Recombinant Human SWSAP1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SWSAP1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SWSAP1(Met1-Pro229) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q6NVH7
Expiration Date	6 months from date of receipt.