

Product Datasheet

Recombinant Human SULT4A1 Protein, No Tag (orb1551586)

Description	Recombinant Human SULT4A1 Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SULT4A1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SULT4A1(Met1-Leu284) was fused without Tag and expressed in E. coli.
Uniprot ID	Q9BR01
Expiration Date	6 months from date of receipt.