

## Product Datasheet

# Recombinant Human SULT1C4 Protein, His Tag (orb1551589)

<b>Description</b>	Recombinant Human SULT1C4 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	SULT1C4
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human SULT1C4(Met1-Lys102) was fused with the N-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="#">O75897</a>
<b>Expiration Date</b>	6 months from date of receipt.