

Product Datasheet

Recombinant Human Stathmin/STMN1 Protein, His Tag (orb1551599)

Description	Recombinant Human Stathmin/STMN1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	STMN1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human STMN1(Ala2-Asp149) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P16949
Expiration Date	6 months from date of receipt.