

Product Datasheet

Recombinant Human SNAP- α /NAPA Protein, His Tag (orb1551627)

Description	Recombinant Human SNAP- α /NAPA Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	NAPA
Note	For research use only
Purity	$\geq 90\%$ as determined by SDS-PAGE
Source	This protein contains the human SNAP- α (Met1-Arg295) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P54920
Expiration Date	6 months from date of receipt.