

Product Datasheet

Recombinant Mouse Serpinf1 Protein, No Tag (orb1551660)

Description	Recombinant Mouse Serpinf1 Protein, No Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Serpinf1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Serpinf1(Asp43-Thr417) was fused without Tag and expressed in E. coli.
Uniprot ID	P97298
Expiration Date	6 months from date of receipt.