

Product Datasheet

Recombinant Human SERPINB6 Protein, His & Trx Tag (orb1551672)

Description	Recombinant Human SERPINB6 Protein, His & Trx Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SERPINB6
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SERPINB6(Met1-Pro376) was fused with the N-terminal His and Trx Tag and expressed in E. coli.
Uniprot ID	P35237
Expiration Date	6 months from date of receipt.