

Product Datasheet

Recombinant Human SelenocysteineLyase/SCLY Protein, His Tag (orb1551704)

Description	Recombinant Human SelenocysteineLyase/SCLY Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SCLY
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SCLY(Met1-Ala445) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q96I15
Expiration Date	6 months from date of receipt.