

Product Datasheet

Recombinant Human SDCBP/Syntenin-1 Protein, His Tag (orb1551707)

Description	Recombinant Human SDCBP/Syntenin-1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SDCBP
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human Sdcbp(Ser2-Val298) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	O00560
Expiration Date	6 months from date of receipt.