

## Product Datasheet

# Recombinant Mouse Sdcbp Protein, His Tag (orb1551708)

<b>Description</b>	Recombinant Mouse Sdcbp Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	Sdcbp
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Sdcbp(Ser2-Val299) was fused with the C-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="#">O08992</a>
<b>Expiration Date</b>	6 months from date of receipt.