

Product Datasheet

Recombinant Human S100A9 Protein, No Tag (orb1551719)

| | |
|------------------------|---|
| Description | Recombinant Human S100A9 Protein, No Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | S100A9 |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human S100A9(Thr2-Pro114) was fused without Tag and expressed in E. coli. |
| Uniprot ID | P06702 |
| Expiration Date | 6 months from date of receipt. |