

Product Datasheet

Recombinant Human S100A8/S100A9 Heterodimer Protein, His Tag (orb1551720)

Description	Recombinant Human S100A8/S100A9 Heterodimer Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	S100A8/S100A9
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human S100A8/S100A9 Heterodimer(Met1-Glu93(S100A8)&Thr2-Pro114(S100A9)) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P06702, P05109
Expiration Date	6 months from date of receipt.