

Product Datasheet

Recombinant Human S100A6 Protein, His Tag (orb1551722)

Description	Recombinant Human S100A6 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	S100A6
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human S100A6(Met1-Gly90) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P06703
Expiration Date	6 months from date of receipt.